

С-Реактивный белок, (CRP-Uhs) Универсальный.

ИттипоTD (Иммунотурбодиметрический с латексными частицами)



Параметры для ввода в программу анализатора iLab-Taurus

Test Name: NN **CRPu** Test Code: IL_ **CRPu-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Type
Serum

Reaction Cycle
☒ Standard ☐ Extended

Reporting Unit
☒ Мг/л* ☐ User Define
Decimal Points: 1
Conv. Factor: 0.000

Methodology
Type
☐ End Point ☒ Rate
Measuring Point
19 - ☒ 29
Photometric
☒ 1 Wavelength ☐ 2 Wavelength
Primary: 510 **Secondary**:
Correction Constant
Slope: 1.000 **Intercept**: 0.000

Test Name: NN **CRPu** Test Code: IL_ **CRPu-d** **Measure** **Sample Reagent** **Ranges Limits** **Calibration**

Sample Volume

	Sample	Dilution	
	Volume	Sample Vol.	Diluent Vol.
1	2.0	0.0	0.0
2	2.0	30.0	60.0
3	4.0	0.0	0.0
4	0.0	0.0	0.0

Diluent:
Diluent Warning Limit: 0 Tests

Reagent Volume

R1: CRPu-d1 **R2**: CRPu-d2

	Volume	Diluent Vol.	stirring	Warning Limit(tests)	Stability(days)
R1	100	0.0	<input checked="" type="checkbox"/>	20	No Control
R2	100	0.0	<input checked="" type="checkbox"/>	20	No Control

Sampling Condition

	Condition No.	1	2	3	4
<input type="checkbox"/>	First Run	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Samp.Vol.Reduction	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Below N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Above N-Range	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Panic H	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	User Range L	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	User Range H	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Noise	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Prozone	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	HIGH!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	ABS!	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>

С-Реактивный белок, (CRP-Uhs) Универсальный.

ИттипоTD (Иммунотурбодиметрический с латексными частицами)



Ranges		Limits																																																																																																								
<table border="1"> <thead> <tr> <th colspan="3">Normal Range</th> </tr> <tr> <th></th> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>Male</td> <td>0.0</td> <td>5.0</td> </tr> <tr> <td>Female</td> <td>0.0</td> <td>5.0</td> </tr> <tr> <td>Other</td> <td>0.0</td> <td>5.0</td> </tr> </tbody> </table>		Normal Range				Lower	Upper	Male	0.0	5.0	Female	0.0	5.0	Other	0.0	5.0	<table border="1"> <thead> <tr> <th colspan="2">Reaction Slope</th> <th colspan="2">Absorbance Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Negative</td> <td><input checked="" type="radio"/> Positive</td> <td><input checked="" type="radio"/> Above</td> <td><input type="radio"/> Below</td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="4"> <table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="4"> <table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </tbody> </table> </td> </tr> </tbody> </table>		Reaction Slope		Absorbance Limit		<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Non-Linear Limit		100	%	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>		Prozone Limit		<input type="radio"/> Above	<input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>			<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	350	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	350	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>				Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0	<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </tbody> </table>				Qualitative		<input type="radio"/> On	<input checked="" type="radio"/> Off
Normal Range																																																																																																										
	Lower	Upper																																																																																																								
Male	0.0	5.0																																																																																																								
Female	0.0	5.0																																																																																																								
Other	0.0	5.0																																																																																																								
Reaction Slope		Absorbance Limit																																																																																																								
<input type="radio"/> Negative	<input checked="" type="radio"/> Positive	<input checked="" type="radio"/> Above	<input type="radio"/> Below																																																																																																							
<table border="1"> <thead> <tr> <th colspan="2">Non-Linear Limit</th> </tr> <tr> <td>100</td> <td>%</td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Non-Linear Limit		100	%	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Prozone Limit</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> Above</td> <td><input checked="" type="radio"/> Below</td> </tr> <tr> <td> <table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td> <table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> </tr> <tr> <td colspan="2"> <table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table> </td> <td colspan="2"></td> </tr> </tbody> </table>		Prozone Limit		<input type="radio"/> Above	<input checked="" type="radio"/> Below	<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>			<table border="1"> <thead> <tr> <th colspan="2">User Range</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		User Range		Lower	Upper	0.01	350	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">ValidRange</th> </tr> <tr> <th>Lower</th> <th>Upper</th> </tr> </thead> <tbody> <tr> <td>0.01</td> <td>350</td> </tr> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		ValidRange		Lower	Upper	0.01	350	<input type="text"/>	<input type="text"/>																																																	
Non-Linear Limit																																																																																																										
100	%																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
Prozone Limit																																																																																																										
<input type="radio"/> Above	<input checked="" type="radio"/> Below																																																																																																									
<table border="1"> <thead> <tr> <th colspan="2">Limit</th> </tr> <tr> <td>0.0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Limit		0.0	<input type="text"/>	<input type="text"/>	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Equation</th> </tr> <tr> <td>none</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>	Equation		none	<input type="text"/>	<input type="text"/>	<input type="text"/>																																																																																													
Limit																																																																																																										
0.0	<input type="text"/>																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
Equation																																																																																																										
none	<input type="text"/>																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
<table border="1"> <thead> <tr> <th colspan="2">Judge Point</th> </tr> <tr> <td>0</td> <td><input type="text"/></td> </tr> </thead> <tbody> <tr> <td><input type="text"/></td> <td><input type="text"/></td> </tr> </tbody> </table>		Judge Point		0	<input type="text"/>	<input type="text"/>	<input type="text"/>																																																																																																			
Judge Point																																																																																																										
0	<input type="text"/>																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
User Range																																																																																																										
Lower	Upper																																																																																																									
0.01	350																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
ValidRange																																																																																																										
Lower	Upper																																																																																																									
0.01	350																																																																																																									
<input type="text"/>	<input type="text"/>																																																																																																									
<table border="1"> <thead> <tr> <th colspan="2">Serum Index Limits</th> </tr> <tr> <th></th> <th>Limit</th> </tr> </thead> <tbody> <tr> <td>Hemolysis</td> <td>0.0</td> </tr> <tr> <td>Icterus</td> <td>0.0</td> </tr> <tr> <td>Lipemia</td> <td>0.0</td> </tr> </tbody> </table>				Serum Index Limits			Limit	Hemolysis	0.0	Icterus	0.0	Lipemia	0.0																																																																																													
Serum Index Limits																																																																																																										
	Limit																																																																																																									
Hemolysis	0.0																																																																																																									
Icterus	0.0																																																																																																									
Lipemia	0.0																																																																																																									
<table border="1"> <thead> <tr> <th colspan="2">Qualitative</th> </tr> </thead> <tbody> <tr> <td><input type="radio"/> On</td> <td><input checked="" type="radio"/> Off</td> </tr> </tbody> </table>				Qualitative		<input type="radio"/> On	<input checked="" type="radio"/> Off																																																																																																			
Qualitative																																																																																																										
<input type="radio"/> On	<input checked="" type="radio"/> Off																																																																																																									

Test Name		Test Code		Measure		Sample Reagent		Ranges Limits		Calibration																																			
NN	CRPu	IL	CRPu-d																																										
<table border="1"> <thead> <tr> <th colspan="2">Calibration</th> </tr> </thead> <tbody> <tr> <td>Method</td> <td>Curve Type</td> </tr> <tr> <td>M-point</td> <td>5P Log</td> </tr> <tr> <td>Repeats</td> <td>Stability</td> </tr> <tr> <td>2</td> <td><input type="text"/></td> </tr> <tr> <td>No</td> <td><input type="text"/></td> </tr> <tr> <td>6</td> <td><input type="text"/></td> </tr> </tbody> </table>				Calibration		Method	Curve Type	M-point	5P Log	Repeats	Stability	2	<input type="text"/>	No	<input type="text"/>	6	<input type="text"/>	<table border="1"> <thead> <tr> <th colspan="2">Condition</th> </tr> </thead> <tbody> <tr> <td><input type="checkbox"/></td> <td>Reagent Blank</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Auto Reagent Blank by Bottle</td> </tr> <tr> <td><input type="checkbox"/></td> <td>Ask for calibration when reagent lot changes</td> </tr> </tbody> </table>				Condition		<input type="checkbox"/>	Reagent Blank	<input type="checkbox"/>	Auto Reagent Blank by Bottle	<input type="checkbox"/>	Ask for calibration when reagent lot changes	<table border="1"> <thead> <tr> <th colspan="2">Calibration/R-Blank Limit</th> </tr> </thead> <tbody> <tr> <td>R-Blank Limit</td> <td>2500.0 mAbs</td> </tr> <tr> <td>Cal Reps Range</td> <td>20.0 %</td> </tr> <tr> <td>Min Cal Reps</td> <td>100.0 mAbs</td> </tr> <tr> <td>Factor Change</td> <td>35 %</td> </tr> <tr> <td>M-Point Curve Fit</td> <td>0.0 %</td> </tr> </tbody> </table>				Calibration/R-Blank Limit		R-Blank Limit	2500.0 mAbs	Cal Reps Range	20.0 %	Min Cal Reps	100.0 mAbs	Factor Change	35 %	M-Point Curve Fit	0.0 %
Calibration																																													
Method	Curve Type																																												
M-point	5P Log																																												
Repeats	Stability																																												
2	<input type="text"/>																																												
No	<input type="text"/>																																												
6	<input type="text"/>																																												
Condition																																													
<input type="checkbox"/>	Reagent Blank																																												
<input type="checkbox"/>	Auto Reagent Blank by Bottle																																												
<input type="checkbox"/>	Ask for calibration when reagent lot changes																																												
Calibration/R-Blank Limit																																													
R-Blank Limit	2500.0 mAbs																																												
Cal Reps Range	20.0 %																																												
Min Cal Reps	100.0 mAbs																																												
Factor Change	35 %																																												
M-Point Curve Fit	0.0 %																																												
<table border="1"> <thead> <tr> <th colspan="3">Calibrator</th> </tr> <tr> <th></th> <th>Calibrator</th> <th>Conc</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>NaCl-0.9</td> <td>0.0</td> </tr> <tr> <td>2</td> <td>TC-CRPu-1</td> <td>**</td> </tr> <tr> <td>3</td> <td>TC-CRPu-2</td> <td>**</td> </tr> <tr> <td>4</td> <td>TC-CRPu-3</td> <td>**</td> </tr> <tr> <td>5</td> <td>TC-CRPu-4</td> <td>**</td> </tr> <tr> <td>6</td> <td>TC-CRPu-5</td> <td>**</td> </tr> </tbody> </table>				Calibrator				Calibrator	Conc	1	NaCl-0.9	0.0	2	TC-CRPu-1	**	3	TC-CRPu-2	**	4	TC-CRPu-3	**	5	TC-CRPu-4	**	6	TC-CRPu-5	**																		
Calibrator																																													
	Calibrator	Conc																																											
1	NaCl-0.9	0.0																																											
2	TC-CRPu-1	**																																											
3	TC-CRPu-2	**																																											
4	TC-CRPu-3	**																																											
5	TC-CRPu-4	**																																											
6	TC-CRPu-5	**																																											

С-Реактивный белок, (CRP-Uhs) Универсальный.

ИттипоTD (Иммунотурбодиметрический с латексными частицами)



* -Рекомендуемый параметр.

** - Вводится из паспорта калибратору.

Калибровка многоточечная (6 точек, 0 + 5 ур.) Spline или 5P Log, калибратор TruCal CRP-U.

В качестве нулевой точки используйте физраствор.

Контроль TruLab Protein уровень 1 и 2 или TruLab CRP уровень 1 и 2.