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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CysC Type: Photometric Version: 1  
 Number: \*\*\* Assay availability: Enabled Date: 18.02.2020  
 Run controls for onboard reagents by: Lot Time: 11:01:08  
 Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: Rate up

Primary Secondary Read times  
 Wavelength: 505 / None Main: 19 -- 28  
 Last required read: 33 Flex: --  
 Absorbance range: -- Color correction: --  
 Sample blank type: None

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CysC Type: Photometric Version: 1  
 Number: \*\*\* Assay availability: Enabled Date: 20.02.2020  
 Run controls for onboard reagents by: Lot Time: 17:14:07  
 Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: CysC  
 Diluent: Saline  
 Diluent dispense mode: Type 0

R1 R2  
 Reagent volume: 135 45  
 Water volume: --  
 Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	1.5				1:1.00	<input checked="" type="radio"/>
1:5	20.0	1.5	80.0		1:5.00	<input type="radio"/>
2x	3.0				2:1.00	<input type="radio"/>

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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CysC Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCalCysC

Calibrator level: Blank: Water Concentration: 0.0000

Cal 1: TruCalCysC1  
Cal 2: TruCalCysC2  
Cal 3: TruCalCysC3  
Cal 4: TruCalCysC4  
Cal 5: TruCalCysC5  
Cal 6: None

Replicates: 2 [Range 1 - 3]

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: CysC Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Spline Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	1.5			
Cal 1: TruCalCysC1	1.5			
Cal 2: TruCalCysC2	1.5			
Cal 3: TruCalCysC3	1.5			
Cal 4: TruCalCysC4	1.5			
Cal 5: TruCalCysC5	1.5			
Cal 6: None				

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