

# ГЛДГ (глютаматдегидрогеназа) (DiaSys)

DGKC-1970



## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type:  Version:   
Number:  Assay availability:  Date:   
Run controls for onboard reagents by:  Time:   
Operator:

☒ Reaction definition ☐ Reagent / Sample ☐ Validity checks

Reaction mode:

Primary Secondary Read times  
Wavelength:  /  Main:  --   
Last required read:  Flex:  --   
Absorbance range:  --  Color correction:  --   
Sample blank type:

Assay of

Done

Cancel

?

Configure assay parameters...

☐ General ☒ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type:  Version:   
Number:  Assay availability:  Date:   
Run controls for onboard reagents by:  Time:   
Operator:

☐ Reaction definition ☒ Reagent / Sample ☐ Validity checks

Reagent:   
Diluent:   
Diluent dispense mode:

R1 R2  
Reagent volume:    
Water volume:   
Dispense mode:

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	: 24				= 1:1.00	<input checked="" type="radio"/>
1:8	: 3				= 1:8.00	<input type="radio"/>
1:20	: 5	16	95		= 1:20.00	<input type="radio"/>

Assay of

Done

Cancel

?

# ГЛДГ (глутаматдегидрогеназа) (DiaSys)

DGKC-1970



## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☐ General ☒ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: GLDH      Assay number: \*\*\*      Date: 20.02.2020  
Calibration method: Linear      Time: 17:14  
Operator: FSE

☒ Calibrators ☐ Volumes ☐ Intervals ☐ Validity checks

Calibrator set: TruCal\_U      Calibrator level: Blank: Water      Concentration: 0.0000

Cal 1: TruCal\_U  
Cal 2: None  
Cal 3: None  
Cal 4: None  
Cal 5: None  
Cal 6: None

Replicates: 2 [Range 1 - 3]

Assay of

Done

Cancel

?

Configure assay parameters...

☐ General ☒ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay: GLDH      Assay number: \*\*\*      Date: 20.02.2020  
Calibration method: Linear      Time: 17:14  
Operator: FSE

☐ Calibrators ☒ Volumes ☐ Intervals ☐ Validity checks

Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	24			
Cal 1: TruCal_U	24			
Cal 2: None				
Cal 3: None				
Cal 4: None				
Cal 5: None				
Cal 6: None				

Assay of

Done

Cancel

?

# ГЛДГ (глутаматдегидрогеназа) (DiaSys)

DGKC-1970



---

Параметры для ввода в программу анализатора Architect c8000 / c4000

Единицы измерения - Ед/л; знаков после запятой 1 (**Configure result units**).

Диапазон метода для сыворотки 2-120 Ед/л (**Configure assay parameters\ Results**).

Референсные интервалы: Муж. 0-7 Ед/л Жен. 0-5 Ед/л (**Configure assay parameters\ Results**).

Калибровка линейная, калибратор TruCal U, нулевой - вода (или физраствор). Возможно использовать калибровку по фактору.

Контроль по TruLab N и TruLab P.