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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Fe Type: Photometric Version: 1  
 Number: 2002 Assay availability: Enabled Date: 18.02.2020  
 Time: 11:01:08  
 Run controls for onboard reagents by: Lot Operator: ADMIN

Reaction definition Reagent / Sample Validity checks

Reaction mode: End up

Primary Secondary Read times  
 Wavelength: 604 / 804 Main: 32 -- 33  
 Last required read: 33  
 Absorbance range: 0.000 -- 2.400 Color correction: --  
 Sample blank type: Self Blank: 14 -- 16

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Fe Type: Photometric Version: 1  
 Number: \*\*\* Assay availability: Enabled Date: 20.02.2020  
 Time: 17:14:07  
 Run controls for onboard reagents by: Lot Operator: FSE

Reaction definition Reagent / Sample Validity checks

Reagent: Fe  
 Diluent: Saline  
 Diluent dispense mode: Type 0

R1 R2  
 Reagent volume: 160 40  
 Water volume: --  
 Dispense mode: Type 2 Type 2

Dilution name	Sample	Diluted sample	Diluent	Water	Dilution factor	Default dilution
Standart	16.0				1:1.00	<input checked="" type="radio"/>
1:2	8.0				1:1.93	<input type="radio"/>
						<input type="radio"/>

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Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Fe Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Linear Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCal\_U  
Calibrator level: Blank: Water Concentration: 0.0000  
Cal 1: TruCal\_U1  
Cal 2: None  
Cal 3: None  
Cal 4: None  
Cal 5: None  
Cal 6: None

Replicates: 2 [Range 1 - 3]

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Fe Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Linear Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

Calibrator level	Sample	Diluted sample	Diluent	Water
Blank: Water	16.0			
Cal 1: TruCal_U1	16.0			
Cal 2: None				
Cal 3: None				
Cal 4: None				
Cal 5: None				
Cal 6: None				

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