


# Лактат (DiaSys)

LDH UV (ферментативный метод с ЛДГ)



## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type: Photometric Version: 1  
Number:  Assay availability:  Date: 18.02.2020  
Run controls for onboard reagents by:  Time: 11:01:08  
Operator: ADMIN

☒ Reaction definition ☐ Reagent / Sample ☐ Validity checks

Reaction mode:

Primary Secondary Read times  
Wavelength:  /  Main:  --   
Last required read:   
Absorbance range:  --  Color correction:  --   
Sample blank type:  Blank:  --

Assay 13 of 54

Done Cancel

?

Configure assay parameters...

☒ General ☐ Calibration ☐ SmartWash ☐ Results ☐ Interpretation

Assay:  Type: Photometric Version: 1  
Number:  Assay availability:  Date: 20.02.2020  
Run controls for onboard reagents by:  Time: 17:14:07  
Operator: FSE

☐ Reaction definition ☒ Reagent / Sample ☐ Validity checks

Reagent:   
Diluent:   
Diluent dispense mode:

R1 R2  
Reagent volume:    
Water volume:    
Dispense mode:

| Dilution name | Sample | Diluted sample | Diluent | Water | Dilution factor | Default dilution                 |
|---------------|--------|----------------|---------|-------|-----------------|----------------------------------|
| Standart      | 2.4    |                |         |       | 1:1.00          | <input checked="" type="radio"/> |
| 1:2           | 20.0   | 4.8            | 80      |       | 1:1.98          | <input type="radio"/>            |
|               |        |                |         |       |                 | <input type="radio"/>            |

Assay of

Done Cancel

?

# Лактат (DiaSys)

LDH UV (ферментативный метод с ЛДГ)



## Параметры для ввода в программу анализатора Architect c8000 / c4000

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Lact Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Linear Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator set: TruCal\_U  
Calibrator level: Blank: Water  
Cal 1: TruCal\_U1  
Cal 2: None  
Cal 3: None  
Cal 4: None  
Cal 5: None  
Cal 6: None  
Concentration: 0.0000  
Replicates: 2 [Range 1 - 3]

Assay of

Done Cancel

?

Configure assay parameters...

General Calibration SmartWash Results Interpretation

Assay: Lact Assay number: \*\*\* Date: 20.02.2020  
Calibration method: Linear Time: 17:14  
Operator: FSE

Calibrators Volumes Intervals Validity checks

Calibrator:

| Calibrator level | Sample | Diluted sample | Diluent | Water |
|------------------|--------|----------------|---------|-------|
| Blank: Water     | 2.4    |                |         |       |
| Cal 1: TruCal_U1 | 2.4    |                |         |       |
| Cal 2: None      |        |                |         |       |
| Cal 3: None      |        |                |         |       |
| Cal 4: None      |        |                |         |       |
| Cal 5: None      |        |                |         |       |
| Cal 6: None      |        |                |         |       |

Assay of

Done Cancel

?

# Лактат (DiaSys)

LDH UV (ферментативный метод с ЛДГ)



---

Параметры для ввода в программу анализатора Architect c8000 / c4000

Единицы измерения - ммоль/л; знаков после запятой 2 (**Configure result units**).

Диапазон метода для сыворотки: 0.1-13.3 ммоль/л (**Configure assay parameters\ Results**).

Референсный интервал: 0.5-2.2 ммоль/л (**Configure assay parameters\ Results**).

При наличии ЛИС, необязательно.

Диапазон нормальных значений указан в соответствии с рекомендациями производителя реагентов.

Калибровка линейная, калибратор TruCal U, нулевой - вода (или физраствор).

Контроль по TruLab N и TruLab P